

Prefilled Syringe Details

The large barrel makes it easy to grip.

(Syringe plunger is pushed in to demonstrate usage)

An audible click signals that the needle guard is activated and is ready for disposal.

The automatic safety device prohibits needle sticks or reuse of the syringe.

One-handed activation of the safety device helps keep hands safely behind the exposed needle.

(Picture pointing out Needle Cover)

(Picture pointing out Needle)

(Picture pointing out Plunger)

(Picture pointing out Safety Device)

RESTART button

[End of transcript]

© 2002-2009 sanofi-aventis U.S. LLC.